

Dee Lawlor

## Introduction to Light Microscopy

Tips and Tricks for Beginners

This book offers a beginner's guide to using light microscopes. It begins with a brief introduction to the physics of optics, which will give the reader a basic grasp of the behaviors of light. In turn, each part of the microscope is explained using clear and simple English, together with detailed photographs and diagrams. The reader will learn the function, care and correct use of each part. A troubleshooting section also helps resolve some of the most common issues encountered in light microscopy.

Most people have a general idea of how to use a microscope, but many never get the full benefit, because they receive no training. With easy-to-follow steps and detailed images, this guide will help everyone achieve the best results, and be confident using their microscope.


This book is intended for anyone using a light microscope, such as university students, people in lab environments, hobbyists, educators who teach science to young children, and anyone with a general interest in these valuable tools.

Life Sciences

ISBN 978-3-030-05392-5



► [springer.com](https://www.springer.com)





<b>1</b>	<b>Introduction</b>	<b>1</b>
1.1	History	2
1.2	The Importance of Glass in Microscopy	6
1.3	What Classifies as a Microscope?	7
1.4	Why We Find Microscopes Appealing	8
	References	9
<b>2</b>	<b>The Advantages of Learning Traditional Microscopy</b>	<b>11</b>
2.1	Skill Acquisition	12
2.2	Enjoying Knowledge	12
2.3	Budget-Friendly	13
2.4	Low Training Requirement	13
2.5	Robust	14
2.6	Reliable	14
2.7	Easy Troubleshooting	15
2.8	The Technology Doesn't Age	15
2.9	Easy to Move	15
	References	16
<b>3</b>	<b>Buying a Microscope</b>	<b>17</b>
3.1	Decide What Is Needed	18
	3.1.1 Who Is Going to Be Using the Microscope?	18
	3.1.2 What Type of Samples Are the Users Working with?	19



3.1.3	How Often Is the Microscope Going to Be Used?	19
3.1.4	Where Is the Microscope Going to Be Used?	19
3.1.5	Will It Stay Set Up or Will It Need to Be Packed Up After Each Use?	20
3.1.6	Does the User Require Image Capture?	20
3.2	Budget and Where to Invest	20
3.3	Microscopes for Children	21
3.4	Buying a Second-Hand Microscope	21
<b>4</b>	<b>The Science of Light</b>	<b>25</b>
4.1	What Is Light?	26
4.2	The Basics of Light	26
4.3	Colours of Light	28
4.4	Movement of Light	29
4.4.1	Refraction	29
4.4.2	Reflection	30
4.4.3	Absorption and Emission	32
4.4.4	Scatter	32
4.5	Photosensitivity: How We Perceive Light	33
4.5.1	The Evolution of Vision	33
4.5.2	Photosensitive Cells	34
4.5.3	Human Vision	34
4.5.4	How We See	35
4.5.5	Colour Vision	35
4.5.6	Variation of Vision Between Individuals	36
4.5.7	Colour-Blindness	36
4.5.8	Variation in Physiology	37
	References	37
<b>5</b>	<b>Introduction to Lenses</b>	<b>39</b>
5.1	Lens Shape and Function	40
5.2	Convex Lenses	41
5.3	Concave Lens	41
5.4	Lens Coating	42
5.5	Anti-reflective Coating	42
5.6	Prisms	43
5.7	Magnification	44
5.8	Numerical Aperture	44



5.9	Resolution	46
5.10	Focus	46
5.11	Aberrations	47
	References	47
<b>6</b>	<b>The Anatomy of the Microscope</b>	<b>49</b>
6.1	Styles of Microscope	50
6.2	The Upright Microscope	50
6.3	The Inverted Microscope	52
6.4	The Stereomicroscope	53
	6.4.1 Types of Stereomicroscope	54
6.5	Petrographic Microscope	54
<b>7</b>	<b>The Eyepiece</b>	<b>57</b>
7.1	Compensating Eyepieces	58
7.2	Eyepoint	58
7.3	Interpupillary Distance	59
7.4	Using Eyepieces as a Glasses Wearer	59
7.5	Diopter Eyepieces	60
	7.5.1 How to Adjust a Diopter Eyepiece	60
7.6	Reticle	61
	7.6.1 How to Calibrate a Reticle	62
7.7	Cleaning an Eyepiece	62
7.8	Scale Bars on Digital Systems	63
	References	63
<b>8</b>	<b>The Objective</b>	<b>65</b>
8.1	How to Read an Objective	67
	8.1.1 Achromatic and Apochromatic	67
	8.1.2 Infinity Correction	68
	8.1.3 Magnification	68
	8.1.4 Numerical Aperture (NA)	68
	8.1.5 Planar	68
	8.1.6 Refractive Index (RI)	70
	8.1.7 Thread Depth and Brand Compatibility	70
	8.1.8 Working Distance (WD)	70
8.2	Air Objective	71
8.3	Oil Immersion Objective	72



8.3.1	How to Use an Air Objective	72
8.3.2	How to Use an Oil Immersion Objective	74
8.4	Cleaning an Objective	77
8.4.1	Cleaning an Air Objective	78
8.4.2	Cleaning an Oil Objective	78
8.5	Summary of Objective Codes	79
	Reference	79
<b>9</b>	<b>The Stage</b>	<b>81</b>
9.1	Focusing the Sample	83
9.2	Slide Navigation Controls	83
9.3	Petrographic Microscope Stage	84
9.4	Stages for Inverted Microscopes	85
<b>10</b>	<b>The Condenser and Diaphragm</b>	<b>87</b>
10.1	The Diaphragm	88
10.1.1	The Disc Diaphragm	89
10.1.2	The Iris Diaphragm	89
10.1.3	How to Centre the Iris Diaphragm	91
10.2	The Condenser	92
10.2.1	How to Focus the Condenser	93
	References	94
<b>11</b>	<b>Light Sources</b>	<b>95</b>
11.1	Different Types of Light Source	96
11.1.1	Internal Light Source	96
11.1.2	External Light Source	97
11.2	Different Types of Bulb	97
11.2.1	LED (Light-Emitting Diodes)	97
11.2.2	Incandescent Bulbs	99
11.2.3	Halogen	100
11.2.4	Mercury	101
11.2.5	Metal Halide	101
11.2.6	Fluorescent	101
11.3	White Balance	102
	References	103



<b>12</b>	<b>Choosing the Right Technique</b>	105
12.1	Brightfield	107
12.1.1	What Does It Look Like?	107
12.1.2	When to Choose Brightfield	107
12.1.3	Tech Specs	108
12.1.4	Sample Requirements	108
12.2	Darkfield	108
12.2.1	What Does It Look Like?	108
12.2.2	When to Choose Darkfield	109
12.2.3	Tech Specs	109
12.2.4	Sample Requirements	110
12.3	Phase Contrast	110
12.3.1	What Does It Look Like?	110
12.3.2	When to Choose Phase Contrast	110
12.3.3	Tech Specs	111
12.3.4	Sample Requirements	111
12.4	Z-Stacking	112
12.4.1	What Does It Look Like?	112
12.4.2	When to Choose Z-Stacking	112
12.4.3	Tech Specs	112
12.4.4	Sample Requirements	113
12.5	Stereomicroscopy	113
12.5.1	What Does It Look Like?	113
12.5.2	When to Choose Stereomicroscopy	113
12.5.3	Tech Specs	113
12.5.4	Sample Requirements	113
<b>13</b>	<b>Sample Prep</b>	115
13.1	Prepared Slide Sets	116
13.2	Equipment for Preparing Samples	116
13.2.1	Glass Slides	116
13.2.2	Coverslips	118
13.2.3	Stains	119
13.2.4	Dry Mount/Wet Mount and Temporary/ Permanent Mount	119
13.3	Health and Safety	120
13.4	Simple Thin-Prep Method for Brightfield Microscopy	120
13.4.1	Cutting the Sample	121
13.4.2	Fixation	122
13.4.3	Staining	123
13.4.4	Coverslipping	123



<b>14</b>	<b>Image Capture</b>	127
14.1	Mounted Camera	128
14.1.1	C-Mounted Cameras	128
14.1.2	Modular Camera	129
14.1.3	Software	130
14.2	Eyepiece Camera	130
14.2.1	Software	131
14.3	Camera Phone Adapter Mount	131
14.3.1	Software	131
14.4	C-Mount Adaptor and a Regular Camera	131
14.5	Imaging Software	132
14.6	Good Practice Points for Image Capture	132
14.7	What Is a Digital Image?	133
14.7.1	How Is a Digital Image Created?	133
14.7.2	What Is a CCD?	133
14.7.3	Pixels	134
14.7.4	Resolution in Digital Images	134
14.7.5	Greyscale	135
14.7.6	RGB Colour	135
14.7.7	Bit Depth	135
14.7.8	Compression	136
14.7.9	Image Formats	137
<b>15</b>	<b>Troubleshooting</b>	139
15.1	A Word of Warning	140
15.2	Problem: There Is No Light Visible in the Eyepieces, or the Light Is Dull or Uneven	141
15.2.1	Power	141
15.2.2	Light Source	141
15.2.3	Diaphragm and Condensers	142
15.2.4	Stage	144
15.2.5	Objective	144
15.2.6	Shutter	145
15.2.7	Eyepieces	146
15.2.8	Camera	146
15.3	Problem: There Is a Shadow in the Image	148
15.4	Problem: Cannot Focus During Air Imaging	148
15.5	Problem: Issues During Oil Imaging	149



15.6	Problem: The Light Is Coming Through the System, But It Is Uneven	150
15.7	Choosing the Right Light Level	150
15.8	Colour of the Sample Is Incorrect	152
15.9	The Wring Test	152
15.10	Centring an Objective on a Petrological Microscope	153
15.11	How to Change the Bulb in an External Light Source	154

<b>Glossary</b>	155
-----------------	-----